


Onsite Labs Aren't Just for Screening Anymore


Timeline for Sample Turnaround



Extraction: 6 minutes using SPME




Analysis: 6 minutes using GC or GC/MS



Reporting: 3 minutes

Onsite Lab Produces Defensible Fixed-Lab-Quality Data



Solid Phase Microextraction (SPME) coupled with fast gas chromatography (GC) delivers sample turnaround times of 15 minutes and daily throughputs of 50 samples per GC, while providing defensible fixed-lab quality.



STONE ENVIRONMENTAL INC

535 Stone Cutters Way
Montpelier, Vermont
05602 USA

Mike Rossi
802.229.2194
www.stone-env.com/profiling